

SQ GUIDE

S Y S T E M



Time
Immediate customized restorations thanks to prosthetically-driven implant planning.



Collaboration
Maximization of treatment outcomes through fully digital collaborative planning



Differentiation
Differentiation of your dental practice with the integrated workflow that enhances your patients' trust and confidence



Efficiency
Efficient surgical procedures thanks to SQ Guided Surgery instruments.


SQ GUIDE

S Y S T E M

Full digital guide system optimized placing SQ Implant


1 Full Solution Guide system

- Place the implant in the correct position by virtually doing the implant placement through the software
- Insert the implant in where you make a plan and have same-day prosthetic setting




2 Exclusive Guide System for SQ Implant

- Computer guided implant placement for SQ Implant that has the advanced SLA surface
- SQ Implant Fixture Size
 - Diameter : $\varnothing 3.5$, $\varnothing 4.0$, $\varnothing 4.5$, $\varnothing 5.0$
 - Length : 7.0mm, 8.0mm, 10.0mm, 12.0mm, 14.0mm



3 SQ GUIDE Sleeve

- Titanium Metal Sleeve
- Sleeve specification
 - Standard offset : 9.0mm
 - Inner diameter : $\varnothing 5.35$
 - Length : 3.5mm, 4.5mm



4 SQ GUIDE Software

- SQ GUIDE is available with DentiqGuide S/W or Implant Studio S/W
- Chair-side production is possible with DENTIS Implant and sleeve library.

