

The New Change 3D Printer  
**ZENITH**



**Head Office (Korea)**

6 Yuram-ro, Dong-gu, Daegu, S. Korea  
Tel +82 53 589 3528 E-mail colin@dentis.co.kr

**DENTIS USA**

6 Centerpointe Drive, Suite 600 La Palma CA 90623  
Tel 323 677 4363 E-mail info@dentisusa.com

## The most beautiful & powerful 3D printer ever

- : A new heated resin tank and sending prints over WiFi bring you the fun of 3d printing.
- : A new user-friendly software applied new UI system enable a user to printer easier than any other machines.

## Professional material specially designe for dentistry

- : Self-developed ZENITH resins enable more reliable, low-cost production of a variety of dental products in-clinic, including surgical guides, dental models, temporary C&B and even castable products.

### Specification



Technology	DLP (Digital Light Processing)
Build Size	128 x 80 x 150 mm
Machine Size (X/Y/Z)	340 x 450 x 430 mm
Layer Thickness	50um, 100um
Light Source	405nm LED
Input file format	Standard STL file
Software	ZENITH Software * Newly introduced
The primary function	<ul style="list-style-type: none"> <li>- 4.3" Touch Screen (Preview Image)</li> <li>- USB Storage</li> <li>- Wifi Connection</li> <li>- Self-Heated Resin Tank</li> <li>- Automatic Firmware Update</li> </ul>

### ZENITH-D Application

: Faster than your 3D printer!



Surgical Guide	50 Minutes / 220 Layers
Dental Mode	40 Minutes / 170 Layers
Temporary C&B	15 Minutes / 120 Layers
Casting Coping	15 Minutes / 120 Layers