


DENTIS SQ | SHORT IMPLANT



Safe and strong fixation

SQ Short Implant




Machined Surface
Strengthening the Bonding force
with soft tissue

Single Thread
Stable depth control

4 Blade Chip Pocket
Maximize Self-Tapping
Capabilities


20°

SQ Short Fixture Line-Up



| Fixture Diameter | Length | Machine Collar | Code No. |
|------------------|--------|----------------|------------|
| 5.0 | 6.0mm | 2.0mm | DSSFR5004S |
| | | 1.0mm | DSSFR5005S |
| | | 0mm | DSSFR5006S |


Fixture Package



Model

Size


Lot no.



Drilling Sequence


Place a fixture 0.5mm to 1mm deep in the lower part of the bone level.

Ø5.0 x 6mm




Normal bone

Hard bone



SQ SHORT&WIDE KIT



DSQS04 DSQS05 DSQS06



DSQSD2206 DSQSTD4006 DSQSTD5006




DSQCSS50 DSQCSW60 DSQCSW70 DSQCSW80

Stopper & Holder

Drill

Countersink



Code - DSQSWK

s-Clean No-Mount Driver & Ratchet Driver

Torque Ratchet

1.25 Hex Driver



DSNDS DRMDSS
DSNDL DRMDSL



DTR



DRHDS125
DRHDL125



Tel. +82-1899-2804 | Fax. +82-53-583-2806
99, Seongseoseo-Ro, Dalseo-Gu, Daegu, Korea



Tel. +1-323-677-4363 | Fax. +1-323-677-4366
6 Centerpointe Drive, Suite 600 La Palma CA 90623



Tel. +021-5111-3828 | Fax. +021-5111-3828
上海市长宁区中山西路933号2205室