

DENTIS CAD/CAM Solution



ScanBody

3shape
exocad


Highly Accurate Library















※ You can download Library at (www.dentis.co.kr) ▶ Product Information ▶ Scanbody.

Link Abutment

3shape
exocad



Diagrams for Digital System

Description	Flow	Tool
Lab Analog	 Digital Lab Analog	
Impression Cap	 Scanbody	
Abutment	 Link  Pre-Milled	 1.25 Hex Driver
	 CoverScrew  Healing Abutment	 1.25 Hex Driver
SQ Fixture	 SQ Ø3.5 ~ 4.0  SQ Ø4.5 ~ 8.0	 No-Mount Driver  Ratchet Driver

ScanBody



Submerged Narrow

Diameter	Length	Type	Code No.
Ø4.0	8mm	Narrow Short	DSSBNSS
	12mm	Narrow Long	DSSBNLS
	16mm	Narrow Extra Long	DSSBNELS

Screw : DSSBSS



Mini

Diameter	Length	Type	Code No.
Ø3.8	8mm	Mini Short	DSSBMSS
	12mm	Mini Long	DSSBMLS
	16mm	Mini Extra Long	DSSBMELS

Screw : DSSBSM



Regular

Diameter	Length	Type	Code No.
Ø4.5	8mm	Regular Short	DSSBRSS
	12mm	Regular Long	DSSBRLS
	16mm	Regular Extra Long	DSSBRELS

Screw : DSSBS



Pre-Milled Abutment

Submerge A-TYPE(Arum)

Narrow

Abutment Diameter	Hex	Code No.
Ø10	Double Hex	DSCBAN10HS
	Non-Hex	DSCBAN10NS

Screw Code : DSASS

Mini

Abutment Diameter	Hex	Code No.
Ø10	Hex 2.1	DSCBAM10HS
	Non-Hex	DSCBAM10NS

Screw Code : 2DSASM

Regular

Abutment Diameter	Hex	Code No.
Ø10	Hex	DSCBA10HS
	Non-Hex	DSCBA10NS
Ø14	Hex	DSCBA14HS
	Non-Hex	DSCBA14NS

Screw Code : DSAS

Screw



Internal A-TYPE(Arum)

Regular

Abutment Diameter	Hex	Code No.
Ø10	Octa 3.1	*DICBA10S
Ø14		*DICBA14S
Ø10	Non-Octa	*DICBA10NS
Ø14		*DICBA14NS

Screw Code : DEAS

* Coming soon



Wide

Abutment Diameter	Hex	Code No.
Ø10	Octa 3.1	*DICBAW10S
Ø14		*DICBAW14S
Ø10	Non-Octa	*DICBAW10NS
Ø14		*DICBAW14NS

Screw Code : DEAS

* Coming soon



Submerge M-TYPE(Manix)

Narrow

Abutment Diameter	Hex	Code No.
Ø10	Double Hex	DSCBMN10HS
	Non-Hex	DSCBMN10NS



Mini

Abutment Diameter	Hex	Code No.
Ø10	Hex 2.1	DSCBMM10HS
	Non-Hex	DSCBMM10NS



Regular

Abutment Diameter	Hex	Code No.
Ø10	Hex	DSCBM10HS
	Non-Hex	DSCBM10NS
Ø14	Hex	DSCBM14HS
	Non-Hex	DSCBM14NS



Internal M-TYPE(Manix)

Regular

Abutment Diameter	Hex	Code No.
Ø10	Octa 3.1	*DICBM10S
Ø14		*DICBM14S
Ø10	Non-Octa	*DICBM10NS
Ø14		*DICBM14NS

* Coming soon



Wide

Abutment Diameter	Hex	Code No.
Ø10	Octa 3.1	*DICBMW10S
Ø14		*DICBMW14S
Ø10	Non-Octa	*DICBMW10NS
Ø14		*DICBMW14NS

* Coming soon



Link Abutment




Link Abutments are used when making Zirconia custom abutments with CAD/CAM equipment

Recommended Torque: 20Ncm (Narrow), 30Ncm (Mini, Regular)

Narrow

Abutment Diameter	Hex	G/H	Code No.
Ø 4.0	HEX	1	DSLAM401045HS
		2	DSLAM402045HS
		3.5	DSLAM403545HS
		5.0	DSLAM405045HS


Screw Code : DSASS



Mini

Abutment Diameter	Hex	G/H	Code No.
Ø 4.0	HEX	1	DSLAM401045HS
		2	DSLAM402045HS
		3.5	DSLAM403545HS
		5.0	DSLAM405045HS


Screw Code : 2DSASM



Regular

Abutment Diameter	Hex	G/H	Code No.
Ø 4.8	HEX	1	DSLAR481045HS
		2	DSLAR482045HS
		3.5	DSLAR483545HS
		5.0	DSLAR485045HS

Screw Code : 2DSAS




Digital Analog

Digital Lab Analogs are used for making models with a 3D printer



Type	Code No.
Narrow	DSDLANS
Mini	DSDLAMS
Regular	DSDLARS



Tighten screw from the bottom

www.dentisimplant.co.kr



Tel. +82-1899-2804 | Fax. +82-53-583-2806
99, Seongseoseo-Ro, Dalseo-Gu, Daegu, Korea



Tel. +1-323-677-4363 | Fax. +1-323-677-4366
6 Centerpointe Drive, Suite 600 La Palma CA 90623



Tel. +021-5111-3828 | Fax. +021-5111-3828
上海市长宁区中山西路933号2205室